

# Universal Smart Box

auditory, visual, & somatosensory  
evoked potentials

& otoacoustic emissions



# our flagship platform

The Universal Smart Box was the first auditory evoked potential and otoacoustic emissions system to use USB connection to connect to a PC. Since then, through a number of design improvements, it has become the most powerful clinical and research platforms for AEP & OAE.


## stellar performance

Over 30 years of engineering design experience, combined with unsurpassed expertise in evoked responses, have culminated in the **next generation bio-amplifier** to bring you **superior data quality** for evoked potentials and otoacoustic emissions.

### Repeatable, reliable data you can count on

- High definition responses
- Cleaner, more robust responses
- Increased signal-to-noise ratio (SNR)
- Lower residual noise

**Reduced test times** without compromising data quality



## the ultimate in flexibility & versatility

The Universal Smart Box can be configured to be an EP system, and OAE system, or an EP/OAE system and can have a maximum of 8 channels. For evoked potentials, the system can be configured to have 1, 2, 4, 6, or 8 channels. OAE systems can have up to two OAE channels. Choose from a variety of add-on modules for the ultimate in flexibility and versatility. Upgrade anytime with minimal or no down time.

## the Smart Box that does it all

### Standard SmartEP modules:

- ECochG
- ABR (click, tone burst, iChirp)
- MLR
- LLR/CAEP

### Standard SmartOAE modules:

- DPOAE
- TEOAE
- SOAE

### Optional SmartEP modules:

- P300/MMN
- eABR
- Chained-Stimuli ABR
- cVEMP, oVEMP
- ASSR
- ENoG
- VEP

## Somatosensory Evoked Potentials

IHS offers an integrated ENoG & Somatosensory EP solution: we manufacture our own Electrical Stimulator Box and Probe.

- Electroneuronography (ENoG) & SSEP on the limbs
- Electrode leads can be connected to surface snap electrodes to deliver stimulus:: no need to hold the probe during the test.
- Low current stimulator option for transtympanic eABR animal research



## Visual Evoked Potentials

Elicit VEP responses using the patterns on a LED light array on the IHS VEP Stimulator.

- Visual patterns can be programmed to stimulate full, half, or quarter fields.
- Stimulator includes centering point to facilitate test subject focusing
- Select from up to 138 different pre-programmed patterns



## for all your research needs

### Advanced options for SmartEP:

- CLAD for high-rate stimulation
- Notched Noise Masking
- Advanced Auditory Research Module
- Complex ABR
- Frequency Following Response
- Acoustic Change Complex
- CHIRP Stimulus Generation Module
- USB Development Kit
- SmartEP-CAM continuous acquisition module
- Auditory P50
- High Frequency ABR up to 32kHz

### Advanced options for SmartOAE:

- Contralateral, ipsilateral, and binaural TEOAE suppression
- Dual OAE probe system
- High Frequency DPOAE for ototoxicity monitoring & animal research





designed for an improved clinical experience

# SmartEP

The ideal clinical tool for recording ECoChG, ABR, and more.

## new and improved user interface simplifies acquisition

- Improved toolbar and button design for fast access to key features
- Easy access to all parameters from a simplified control panel and streamlined menus
- Quickly load your own or preset protocols
- Easily view ongoing EEG display for quick assessment of patient state during testing
- Great variety of options allow you to perform the tests the way you want
- Choose from a variety of stimuli, or generate or import your own custom stimuli
- Display or hide a subtle vertical grid or horizontal baseline



### streamlined workflow


- Easily mark waveforms using over thirty pre-defined peak labels, or create your own custom labels. Easily adjust them using a mouse or keyboard
- View latencies and amplitudes of peaks directly on waveforms and in newly embedded recording information panel.
- Automatically arrange recordings by intensity, acquisition order, stimulus frequency, or rate
- Quickly resize the waveforms using the zoom in/out buttons
- Multi-page display and reports
- Easy PDF report generation
- Auto-save reports on program exit

### iChirp™ stimuli included

The intelligent Chirp for SmartEP and SmartEP-ASSR is included in the base package.

- Broadband and narrowband (500, 1000, 2000, 4000 Hz)
- Improved threshold detection
- Robust amplitude responses
- Optimized wave V identification
- Optional, innovative custom chirp design utility


### Beneficial for recording ABR in awake and active patients.



2000 Hz tone burst ABR (left) vs 2000 Hz iChirp ABR (right)



### smart features

- Change most test parameters with a single click
- Set your own display scale
- Easily mark waveforms using over thirty pre-defined peak labels, or create your own custom labels
- View latency and amplitude of peaks directly on the waveform
- Latency-Intensity graphs indicating normative data ranges are automatically generated from marked waveforms
- Quickly add, subtract, invert, time shift, or cross-correlate recordings
- Cross-correlation value displayed in information panel




### superior electrocochleography

Our next generation amplifiers combined with the non-invasive IHS Lilly TM-Wick Electrodes produce more robust and repeatable ECoChGs.




- View repeatability, SNR, and Residual Noise value for each recording
- Split-sweep view to visualize single recording repeatability
- Export to ASCII single recording or page



Improved SP/AP amplitude and area curve ratio analysis and automatic calculation.

### smarter averaging display options

- Option to automatically acquire and store data in sweep blocks for more powerful processing
- Easily analyze acquired waveforms using additional averaging techniques for further noise reduction
- Averaging techniques includes traditional linear, median, and weighted



- Objective response measurements provide indicators of recording quality
- Automated averaging stopping rules using residual noise measures allow data to be acquired automatically with consistent quality and noise levels



# SmartDPOAE

Screening and diagnostic distortion product otoacoustic emissions.


- Fast and easy setup with up to 41 frequencies per ear in a single test
- Automatic probe-fit check and in-ear calibration for increased accuracy
- Easy-to-interpret colorful DP Grams and detailed information for each frequency tested
- Clear Pass or Refer indications based on user-selected passing criteria
- User-customizable display of normative ranges on the DP Gram facilitates response analysis
- High frequency option for ototoxicity monitoring
- Built-in scripting feature allows you to define sequences of frequencies and intensities for automated data collection
- Optional graphical display of noise standard deviation for improved interpretation



# SmartEP-ASSR

Full-featured screening and diagnostic Auditory Steady State Response System.


- Provides quick, accurate threshold detection using automated statistical analysis
- Test both ears at the same time, four frequencies per ear
- iChirp (broadband & frequency specific) for robust amplitudes and harmonic component analysis for improved threshold detection and reduced test times
- Automated audiogram generation in SPL and HL
- Cost effective add-on to SmartEP



# SmartTrOAE

Screening and diagnostic transient evoked and spontaneous otoacoustic emissions.

- Fast and easy test setup and data analysis
- Automatic probe-fit check and in-ear calibration
- Clear Pass or Refer indications based on user-selected passing criteria
- Displays of the OAE time signal, frequency analysis and the ear canal response
- Use clicks, tones, or user-defined stimulus files
- Time-Frequency plots can be used to illustrate how the frequency composition of transient OAE responses, Noise, and SNR change over time



Available in the Universal Smart Box is a dual-probe option that allows for the acquisition of contralateral, ipsilateral, and binaural TEOAE suppression recordings. This option includes a Suppression Analysis module for temporal and spectral comparison of control and suppression data.

# Intelligent VRA

Automated visual reinforcement audiometry.

- Increased reliability & accuracy by a single examiner
- Choose from our variety of 4 and 10 second colorful, animated wide-screen video clips, or use your own video clips
- Use one of our three automated test routines, administer a speech discrimination paradigm (IVRISD), or run a VRA test manually
- Probe Trials maintain patient attention while testing near threshold
- Control Trials allow you to determine the reliability of a test
- Trial-by-trial reports include detailed information for each test sequence
- Final report includes audiogram and threshold for each frequency tested



**CAST™**  
Classification of Audiograms by Sequential Testing selects the best-fitting audiogram from 9 patterns, for fast and efficient screening.

**OHTA™**  
Optimized Hearing Test Algorithm is designed to test four frequencies, non-sequentially in an intensity staircase fashion.

**5-up/5-down**  
Automated 'step-up, step-down' intensity staircase procedure for testing thresholds at up to four selected frequencies.

# Specifications

## SmartEP

Adjustable Gain: 5K - 300K  
Adjustable High Pass and Low Pass filters (-6 dB/Oct)  
Stimulus: Clicks, Tones, iChirps, Complex, and user-defined files  
Stimulus duration in  $\mu$ sec or cycles  
Stimulus Envelopes: Rectangular, Blackman, Cosine, Hamming, Hann, Bartlett, Trapezoidal (Rise/fall time), Extended Cosine (Rise/fall time), Triangular, Gaussian  
Stimulus presented continuously or only while acquiring  
Ipsilateral and Contralateral noise masking. Specified level or tracking the stimulus level

## SmartEP-ASSR

Gain: 100K  
High Pass Filter: 30 Hz  
Low Pass Filter: 300 Hz  
Stimulus: Clicks, Tones, iChirps, and user-defined files  
Frequencies: 250, 500, 1000, 2000, 4000, and 8000 Hz  
Simultaneous testing of both ears  
Test up to four frequencies per ear

## SmartDPOAE

Up to 41 frequencies per DPGram  
DP I/O Function

## SmartTrOAE

Response window: 300 - 6000 Hz  
Stimulus: Clicks and Tones  
Contralateral, Ipsilateral, and Binaural suppression  
Dual probe option

## EP Amplifier

Two channels  
A/D Converter: 16-bit  
Sampling rate: 200 to 40000 Hz  
High Pass: 1 - 500 Hz  
Low Pass: 30 - 5000 Hz  
Adjustable artifact rejection level and time region  
Line Frequency Notch Filter (-12 dB/Oct)  
Common Mode Rejection:  
 $\geq 105$  dB @ 1 kHz  
 $\geq 120$  dB @ 60/50 Hz, notch filter off  
Noise Level:  $\leq 0.333$   $\mu$ V RMS  
Input Impedance:  $> 5$  MOhms

## Transducers

ER-3C Insert Earphones:  
Intensity: 0 - 130 dB SPL  
Frequency Range: 125 - 10000 Hz  
ER-2 Insert Earphones:  
Intensity: 0 - 118 dB SPL  
Frequency Range: 125 - 16000 Hz  
High Frequency Transducers:  
Intensity: 0 - 94 dB SPL  
Frequency Range: 2000 - 32000 Hz  
TDH Headphones:  
Intensity: 0 - 120 dB SPL  
Frequency Range: 125 - 12000 Hz  
Bone Conductor:  
Intensity: 0 - 98 dB SPL  
Frequency Range: 250 - 8000 Hz  
ER-10D OAE Probe:  
Intensity: 0 - 100 dB SPL  
Frequency Range: 125 - 16000 Hz  
ER-10B+ OAE Microphone:  
Sound field amplifier and speakers  
Auxiliary output channel for ipsilateral masking and stimulus mixing

## Power Requirements

115 - 230 VAC, 50/60 Hz,  
560 - 350 mA, 30 W

## Operating Environment

Portable Equipment  
Indoor use  
Operating temperature: 15 °C - 40 °C  
Relative humidity: 15% to 90% at 40 °C non-condensing  
Altitude: 0 - 3000 m

## Storage

Temperature: 0 °C - 50 °C  
Atmospheric Pressure: none specified

## Standards Compliance

Safety: IEC 60601-1 Class II, Type BF  
EMC: IEC 60601-1-2  
EP: IEC 60601-2-40  
Medical Device Directive: 93/42/EEC

## Computer Requirements

Windows 10 operating system  
Minimum 4 GB RAM  
Minimum 5 GB hard drive space  
Min display vertical resolution of 900 px, Full HD recommended.  
Grounded, 3-prong power supply  
Compliant with IEC 60950  
Mouse or other pointing device  
One available USB Port  
Removable media, network drive, or secure Internet storage site for data backup (recommended)  
Printer (optional)

\*May not be available in all markets